



Product Specification
**Screw cap tube, 30 ml, (LxØ): 107 x 25 mm, PP,
with print**



Product description

Order number	62.543
Product description	Screw cap tube, working volume: 30 ml, (LxØ): 107 x 25 mm, material: PP, skirted conical base, transparent, screw cap, natural, cap assembled, with print, label/print: black, with graduations, 125 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	black
Graduation	yes
Closure type	screw cap
Base shape	skirted conical base
Volume of work	30 ml
Centrifugation max (RCF)	8800 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 30 ml, (LxØ): 107 x 25 mm, PP,
with print**

Page 2

Size

Diameter	25 mm
Length of product	107 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

Purity & certification

Satisfies the requirement	ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	125
Piece (s) / outer case	500
Piece(s) / pallet	12000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0122 kg
Weight of case	6.84 kg
EAN of inner box	4038917082365
EAN case	4038917001465

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT