



### Product description

Order number	62.547.004
Product description	Screw cap tube, working volume: 50 ml, (LxØ): 114 x 28 mm, material: PP, conical base, transparent, HD-PE, cap: red, cap assembled, with print, label/print: blue/white, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 25 piece(s)/polystyrene container

### Product characteristics

Label/ Print	with print
Colour of print/label	blue/white
Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Base shape	conical base
Volume of work	50 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP,  
with print**

**Size**

Diameter 28 mm  
Length of product 114 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap red

**Purity & certification**

Satisfies the requirement ADR  
Product category In vitro diagnostic device, CE  
CE certified CE  
Purity standard free from DNA/DNase/RNase, free from  
pyrogens/endotoxins, non-cytotoxic, sterile  
Sterilisation Electron irradiation  
Pyrogen-free yes  
Batched yes

**Packaging**

Minimum order qty. 300  
Type of smallest polystyrene container  
subpackaging  
Piece(s) / inner box 25  
Piece (s) / outer case 300  
Piece(s) / pallet 8400  
Depth of case 588 mm  
Width of case 398 mm  
Height of case 255 mm  
Case volume 0.0597 cbm  
Weight of product 0.0128 kg  
Weight of case 4.83 kg  
EAN of inner box 4038917111829  
EAN case 4038917013413

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

