



Product Specification
**Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP,
with print**



Product description

Order number	62.547.274
Product description	Screw cap tube, working volume: 50 ml, (LxØ): 114 x 28 mm, material: PP, conical base, transparent, HD-PE, cap: natural, cap assembled, with print, label/print: blue/white, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 25 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	blue/white
Graduation	yes
Closure type	screw cap
Base shape	conical base
Volume of work	50 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP,
with print**

Page 2

Size

Diameter 28 mm
Length of product 114 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap natural

Purity & certification

Satisfies the requirement ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard free from DNA/DNase/RNase, free from
pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation Electron irradiation
Pyrogen-free yes
Batched yes

Packaging

Minimum order qty. 300
Type of smallest bag
subpackaging
Piece(s) / inner box 25
Piece (s) / outer case 300
Piece(s) / pallet 9000
Depth of case 580 mm
Width of case 330 mm
Height of case 255 mm
Case volume 0.0488 cbm
Weight of product 0.0128 kg
Weight of case 4.61 kg
EAN of inner box 4038917242561
EAN case 4038917242554

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT