




**Product Specification**  
**Screw cap tube, 10 ml, (LxØ): 79 x 16 mm, PP,  
with print**

**Page 1**



### Product description

Order number	62.551.011
Product description	Screw cap tube, working volume: 10 ml, (LxØ): 79 x 16 mm, material: PP, round base with skirted base, transparent, screw cap, natural, cap enclosed, with print, label/print: white, with graduations, 500 piece(s)/bag

### Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	screw cap
Base shape	round base with skirted base
Volume of work	10 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 10 ml, (LxØ): 79 x 16 mm, PP,  
with print**

**Page 2**

**Size**

Diameter 16 mm  
Length of product 79 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap natural

**Purity & certification**

Satisfies the requirement IATA  
Product category no medical device | IVD  
Batched yes

**Packaging**

Minimum order qty. 1000  
Type of smallest subpackaging bag  
Piece(s) / inner box 500  
Piece (s) / outer case 1000  
Piece(s) / pallet 30000  
Case volume 0.046 cbm  
Weight of product 0.0048 kg  
Weight of case 5.35 kg  
EAN of inner box 4038917112772  
EAN case 4038917033404

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**