




Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS,
with print**



Product description

Order number	62.553.002
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PS, conical base, transparent, screw cap, red, cap assembled, with print, label/print: white/blue, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 50 piece(s)/rack

Product characteristics

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Closure type	screw cap
Base shape	conical base
Volume of work	15 ml
Centrifugation max (RCF)	4100 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS,
with print**

Size

Diameter 17 mm
Length of product 120 mm

Material & colours

Product material Polystyrene (PS)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap red

Purity & certification

Satisfies the requirement IATA
Product category no medical device | IVD
Purity standard free from DNA/DNase/RNase, free from
pyrogens/endotoxins, non-cytotoxic, sterile
Non-cytotoxic yes
Batched yes

Packaging

Minimum order qty. 500
Type of smallest rack
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 18000
Case volume 0.0436 cbm
Weight of product 0.0101 kg
Weight of case 5.31 kg
EAN of inner box 4038917112949
EAN case 4038917033435

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

