



**Product Specification**  
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS,  
with print**



**Product description**

Order number	62.553.041
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PS, conical base, transparent, screw cap, yellow, cap assembled, with print, label/print: white/blue, with graduations, 50 piece(s)/bag

**Product characteristics**

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Closure type	screw cap
Base shape	conical base
Volume of work	15 ml
Centrifugation max (RCF)	4100 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS,  
with print**

**Size**

Diameter	17 mm
Length of product	120 mm

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow

**Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	20000
Depth of case	428 mm
Width of case	388 mm
Height of case	225 mm
Case volume	0.0374 cbm
Weight of product	0.0077 kg
Weight of case	4.35 kg
EAN of inner box	4038917112963
EAN case	4038917013857

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

