




Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS,
with print**



Product description

| | |
|---------------------|---|
| Order number | 62.553.205 |
| Product description | Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PS, conical base, transparent, screw cap, red, cap assembled, with print, label/print: white/blue, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 50 piece(s)/bag |

Product characteristics

| | |
|--------------------------|--------------|
| Label/ Print | with print |
| Colour of print/label | white/blue |
| Graduation | yes |
| Closure type | screw cap |
| Base shape | conical base |
| Volume of work | 15 ml |
| Centrifugation max (RCF) | 4100 x g |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS,
with print**

Size

Diameter 17 mm
Length of product 120 mm

Material & colours

Product material Polystyrene (PS)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap red

Purity & certification

Satisfies the requirement IATA
Product category no medical device | IVD
Purity standard free from DNA/DNase/RNase, free from
pyrogens/endotoxins, non-cytotoxic, sterile
Non-cytotoxic yes
Batched yes

Packaging

Minimum order qty. 500
Type of smallest bag
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 24000
Case volume 0.0366 cbm
Weight of product 0.0084 kg
Weight of case 4.45 kg
EAN of inner box 4038917112987
EAN case 4038917033442

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

