



Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP,
with print**



Product description

Order number	62.554.002
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PP, conical base, transparent, screw cap, red, cap assembled, with print, label/print: white/blue, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 50 piece(s)/polystyrene container

Product characteristics

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Base shape	conical base
Volume of work	15 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP,
with print**

Size

Diameter 17 mm
Length of product 120 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap red

Purity & certification

Satisfies the requirement IATA
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard free from DNA/DNase/RNase, free from
pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation Electron irradiation
Pyrogen-free yes
Batched yes

Packaging

Minimum order qty. 500
Type of smallest polystyrene container
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 16000
Depth of case 435 mm
Width of case 143 mm
Height of case 702 mm
Case volume 0.0436 cbm
Weight of product 0.0057 kg
Weight of case 4.13 kg
EAN of inner box 4038917565103
EAN case 4038917565110

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

