



Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP,
with print**

Page 1



Product description

Order number	62.554.205
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PP, conical base, transparent, screw cap, red, cap assembled, with print, label/print: white/blue, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 50 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Base shape	conical base
Volume of work	15 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP,
with print**

Size

Diameter	17 mm
Length of product	120 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Purity standard	free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation	Electron irradiation
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	24000
Depth of box	300 mm
Width of box	340 mm
Height of box	100 mm
Depth of case	429 mm
Width of case	391 mm
Height of case	219 mm
Case volume	0.0366 cbm
Weight of product	0.0073 kg
Weight of case	4 kg
EAN of inner box	4038917113144
EAN case	4038917033497

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

