



**Product Specification**  
**Screw cap tube, 50 ml, (LxØ): 115 x 28 mm, PP,  
with print**



**Product description**

|                     |   |
|---------------------|---|
| Order number        | 62.559.205  |
| Product description | Screw cap tube, working volume: 50 ml, (LxØ): 115 x 28 mm, material: PP, skirted conical base, transparent, screw cap, red, cap assembled, with print, label/print: white/blue, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 25 piece(s)/bag |

**Product characteristics**

|                          |                      |
|--------------------------|----------------------|
| Label/ Print             | with print           |
| Colour of print/label    | white/blue           |
| Graduation               | yes                  |
| Closure type             | screw cap            |
| Base shape               | skirted conical base |
| Volume of work           | 50 ml                |
| Centrifugation max (RCF) | 8000 x g             |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 50 ml, (LxØ): 115 x 28 mm, PP,  
with print**

**Size**

Diameter 28 mm  
Length of product 115 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap red

**Purity & certification**

Satisfies the requirement IATA  
Product category no medical device | IVD  
Purity standard free from DNA/DNase/RNase, free from  
pyrogens/endotoxins, non-cytotoxic, sterile  
Sterilisation Electron irradiation  
Pyrogen-free yes  
Non-cytotoxic yes  
Batched yes

**Packaging**

Minimum order qty. 300  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 25  
Piece (s) / outer case 300  
Piece(s) / pallet 8400  
Case volume 0.0649 cbm  
Weight of product 0.0155 kg  
Weight of case 5.19 kg  
EAN of inner box 4038917376402  
EAN case 4038917376396

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

