



## Product Specification

**Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm,  
conical false bottom, rounded tube bottom, PP,  
cap assembled, 100 piece(s)/bag**

**Page 1**



### Product description

Order number	62.612
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 92 x 15.3 mm, conical false bottom, rounded tube bottom, transparent, material: PP, with print, label/print: white, with graduations, cap assembled, 100 piece(s)/bag, 1,000 piece(s)/carton

### Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Base shape	conical false bottom, rounded tube bottom
Volume of work	5 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-09-01

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm,  
conical false bottom, rounded tube bottom, PP,  
cap assembled, 100 piece(s)/bag**

**Page 2**

### Size

Diameter	15.3 mm
Length of product	92 mm
Length excluding cap	92 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.005 kg
Weight of case	5.69 kg
EAN of inner box	4038917082884
EAN case	4038917001908

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-09-01

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**