



**Product Specification**  
**Screw cap tube, 2.5 ml, (LxØ): 75 x 13 mm,  
conical false bottom, rounded tube bottom, PP,  
425 piece(s)/StackPack**



**Product description**

Order number	63.1575
Product description	Screw cap tube, working volume: 2.5 ml, (LxØ): 75 x 13 mm, conical false bottom, rounded tube bottom, transparent, material: PP, with paper thermo label, 425 piece(s)/StackPack, 3,400 piece(s)/carton

**Product characteristics**

Label/ Print	with paper thermo label
Closure type	screw cap
Base shape	conical false bottom, rounded tube bottom
Volume of work	2,5 ml

**Size**

Diameter	13 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 2.5 ml, (LxØ): 75 x 13 mm,  
conical false bottom, rounded tube bottom, PP,  
425 piece(s)/StackPack**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	3400
Type of smallest subpackaging	StackPack
Piece(s) / inner box	425
Piece (s) / outer case	3400
Piece(s) / pallet	102000
Depth of box	381 mm
Width of box	77 mm
Height of box	176 mm
Depth of case	398 mm
Width of case	358 mm
Height of case	325 mm
Case volume	0.0463 cbm
Weight of product	0.0026 kg
Weight of case	9.91 kg
EAN of inner box	4038917313506
EAN case	4038917313490

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**