



Product Specification
**Screw cap tube, 5 ml, 5 ml, (LxØ): 92 x 15.3 mm,
conical false bottom, rounded tube bottom, PP,
250 piece(s)/StackPack**

Product description

Order number 63.1592
Product description Screw cap tube, working volume: 5 ml, nominal volume: 5 ml, (LxØ): 92 x 15.3 mm, conical false bottom, rounded tube bottom, transparent, material: PP, with paper thermo label, 250 piece(s)/StackPack, 1,000 piece(s)/carton

Product characteristics

Label/ Print with paper thermo label
Closure type screw cap
Base shape conical false bottom, rounded tube bottom
Volume of work 5 ml

Size

Sample volume 5 ml
Diameter 15.3 mm
Length of product 92 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Tube material Polypropylene (PP)

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest StackPack
subpackaging
Piece(s) / inner box 250
Piece (s) / outer case 1000
Piece(s) / pallet 48000
Depth of box 336 mm
Width of box 95 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**Screw cap tube, 5 ml, 5 ml, (LxØ): 92 x 15.3 mm,
conical false bottom, rounded tube bottom, PP,
250 piece(s)/StackPack**

Page 2

Height of box	177 mm
Depth of case	358 mm
Width of case	348 mm
Height of case	207 mm
Case volume	0.0258 cbm
Weight of product	0.0042 kg
Weight of case	4.93 kg
EAN of inner box	4038917269452
EAN case	4038917269445

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT