



**Product Specification**  
**Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm,  
conical false bottom, rounded tube bottom, PP,  
250 piece(s)/StackPack**

**Product description**

Order number	63.1593
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 92 x 15.3 mm, conical false bottom, rounded tube bottom, transparent, material: PP, with paper thermo label, 250 piece(s)/StackPack, 1,000 piece(s)/carton

**Product characteristics**

Label/ Print	with paper thermo label
Closure type	screw cap
Base shape	conical false bottom, rounded tube bottom
Volume of work	5 ml

**Size**

Diameter	15.3 mm
Length of product	92 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	StackPack
Piece(s) / inner box	250
Piece (s) / outer case	1000
Piece(s) / pallet	48000
Depth of box	336 mm
Width of box	86 mm
Height of box	171 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm,**  
**conical false bottom, rounded tube bottom, PP,**  
**250 piece(s)/StackPack**

**Page 2**

Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0036 kg
Weight of case	4.23 kg
EAN of inner box	4038917087179
EAN case	4038917005647

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**