



Product Specification
**Screw cap tube, 5 ml, (LxØ): 57 x 16.5 mm, PP,
with paper label**



Product description

Order number	63.542.709
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 57 x 16.5 mm, material: PP, flat base, transparent, screw cap, natural, cap enclosed, with paper label, label/print: white/black, with graduations, 1,000 piece(s)/bag

Product characteristics

Label/ Print	with paper label
Colour of print/label	white/black
Graduation	yes
Closure type	screw cap
Base shape	flat base
Volume of work	5 ml
Centrifugation max (RCF)	4000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap tube, 5 ml, (LxØ): 57 x 16.5 mm, PP,
with paper label

Page 2

Size

Diameter	16.5 mm
Length of product	57 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

Purity & certification

Satisfies the requirement	IATA
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	42000
Depth of case	388 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0367 cbm
Weight of product	0.0043 kg
Weight of case	4.73 kg
EAN of inner box	4038917111621
EAN case	4038917013277

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT