



**Product Specification**  
**Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, PP,**  
**with paper label**



**Product description**

Order number	63.9922.252
Product description	Screw cap tube, working volume: 25 ml, (LxØ): 90 x 25 mm, material: PP, skirted conical base, transparent, screw cap, white, cap assembled, with paper label, label/print: white/black, Patient information, with graduations, 500 piece(s)/bag

**Product characteristics**

Label/ Print	with paper label
Colour of print/label	white/black
Graduation	yes
Closure type	screw cap
Base shape	skirted conical base
Volume of work	25 ml
Centrifugation max (RCF)	4500 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, PP,**  
**with paper label**

**Page 2**

**Size**

Diameter	25 mm
Length of product	90 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white

**Purity & certification**

Satisfies the requirement	IATA
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	500
Piece(s) / pallet	15000
Depth of case	608 mm
Width of case	324 mm
Height of case	264 mm
Case volume	0.052 cbm
Weight of product	0.0084 kg
Weight of case	4.98 kg
EAN of inner box	4038917151559
EAN case	4038917027779

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**