



Product Specification

**Multi-purpose container, max. volume: 120 ml,
(LxØ): 105 x 44 mm, graduated, PP, transparent**

Page 1



Product description

Order number	75.9922.421
Product description	Multi-purpose container, max. volume: 120 ml, (LxØ): 105 x 44 mm, transparent, cap: transparent, graduated, material: PP, screw cap, cap assembled, 250 piece(s)/bag

Product characteristics

Graduation	yes
Closure type	screw cap
Base shape	flat base
Volume of work	120 ml

Size

Sample volume	120 ml
Max. volume	120 ml
Diameter	44 mm
Height of product	105 mm
Length of product	105 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Multi-purpose container, max. volume: 120 ml,
(LxØ): 105 x 44 mm, graduated, PP, transparent**

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	250
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	250
Piece(s) / pallet	4000
Depth of case	548 mm
Width of case	408 mm
Height of case	373 mm
Case volume	0.0834 cbm
Weight of product	0.0188 kg
Weight of case	5.62 kg
EAN of inner box	4038917151870
EAN case	4038917038478

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT