




## Product Specification

**Multi-purpose container, max. volume: 250 ml,  
(LxØ): 78 x 70 mm, graduated, PP, translucent,  
with paper label**

**Page 1**



### Product description

Order number	75.9922.533
Product description	Multi-purpose container, max. volume: 250 ml, (LxØ): 78 x 70 mm, translucent, cap: natural, graduated, material: PP, screw cap, cap assembled, with paper label, label/print: black, 80 piece(s)/bag

### Product characteristics

Label/ Print	with paper label
Colour of print/label	black
Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Base shape	flat base
Volume of work	250 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Multi-purpose container, max. volume: 250 ml,  
(LxØ): 78 x 70 mm, graduated, PP, translucent,  
with paper label**

**Page 2**

### Size

Sample volume	250 ml
Max. volume	250 ml
Diameter	70 mm
Height of product	78 mm
Length of product	78 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	translucent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural
Tube material	Polypropylene (PP)

### Purity & certification

Satisfies the requirement	sealed at 9"hg/10min
Latex-free	yes
Pyrogen-free	no
Safety data sheet	no

### Packaging

Minimum order qty.	240
Type of smallest subpackaging	bag
Piece(s) / inner box	80
Piece (s) / outer case	240
Piece(s) / pallet	3840
Depth of case	646 mm
Width of case	498 mm
Height of case	450 mm
Case volume	0.1445 cbm
Weight of case	10.93 kg
EAN of inner box	4038917152112
EAN case	4038917235013

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**