



### **Product description**

Order number	65.716.303
Product description	Screw cap, red, suitable for screw cap micro tubes, 500 piece(s)/bag

### **Product characteristics**

Cap	flat cap
Closure type	screw cap
Suitable for	screw cap micro tubes

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	red
O-ring material	Ethylene propylene diene rubber (EPDM)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap, red, suitable for screw cap micro tubes**

**Purity & certification**

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	10000
Piece(s) / pallet	720000
EAN of inner box	4038917925204
EAN case	4038917925211

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

