



Product description

Order number	72.690.375
Product description	Reaction tube, working volume: 1.5 ml, material: PP, blue, cap: blue, retaining loop cap, cap attached, with injected graduation and writing space, 500 piece(s)/bag

Product characteristics

Graduation	yes
Cap	retaining loop cap
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g
Autoclavable (max. °C)	121 °C

Size

Sample volume	1.5 ml
Length of product	39 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	blue
Colour of cap	blue

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	160000
Case volume	0.0558 cbm
Weight of product	0.001 kg
Weight of case	5.6 kg
EAN of inner box	4038917369442
EAN case	4038917369435

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

