



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 72.691.007   |
| Product description | Reaction tube, working volume: 2 ml, material: PP, violet, cap: violet, retaining loop cap, cap attached, with injected graduation and writing space, 500 piece(s)/bag |

### **Product characteristics**

|                          |                    |
|--------------------------|--------------------|
| Graduation               | yes                |
| Cap                      | retaining loop cap |
| Closure type             | hinged lid         |
| Volume of work           | 2 ml               |
| Centrifugation max (RCF) | 20000 x g          |

### **Size**

|                   |         |
|-------------------|---------|
| Diameter          | 10.8 mm |
| Length of product | 40 mm   |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | violet             |
| Closure material  | Polypropylene (PP) |
| Colour of cap     | violet             |

### **Purity & certification**

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Batched          | yes                           |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 5000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 5000          |
| Piece(s) / pallet             | 120000        |
| Depth of case                 | 398 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 415 mm        |
| Case volume                   | 0.0657 cbm    |
| Weight of product             | 0.0011 kg     |
| Weight of case                | 6.06 kg       |
| EAN of inner box              | 4038917312974 |
| EAN case                      | 4038917312967 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

