



### Product description

Order number	72.701.201
Product description	SafeSeal reaction tube, working volume: 5 ml, material: PP, transparent, cap: natural, SafeSeal cap, cap attached, with injected graduation and writing space, Biosphere® plus, 1 piece(s)/blister

### Product characteristics

Graduation	yes
Cap	SafeSeal cap
Volume of work	5 ml
Centrifugation max (RCF)	25000 x g

### Size

Diameter	17 mm
Length of product	60 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Batched	yes

### **Packaging**

Minimum order qty.	50
Type of smallest subpackaging	blister
Piece(s) / inner box	25
Piece (s) / outer case	50
Piece(s) / pallet	7200
Depth of box	201 mm
Width of box	151 mm
Height of box	126 mm
Depth of case	323 mm
Width of case	213 mm
Height of case	145 mm
Case volume	0.01 cbm
Weight of product	0.0031 kg
Weight of case	0.51 kg
EAN of inner box	4038917920506
EAN case	4038917920513

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

