



### **Product description**

Order number	72.701.600
Product description	SafeSeal reaction tube, working volume: 5 ml, material: PP, transparent, cap: natural, SafeSeal cap, cap attached, with injected graduation and writing space, PCR Performance Tested, Low protein-binding, 25 piece(s)/bag

### **Product characteristics**

Graduation	yes
Cap	SafeSeal cap
Volume of work	5 ml
Centrifugation max (RCF)	25000 x g

### **Size**

Diameter	17 mm
Length of product	60 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

### **Packaging**

Minimum order qty.	200
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	200
Piece(s) / pallet	15600
Depth of box	189 mm
Width of box	151 mm
Height of box	97 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.003 kg
Weight of case	1.18 kg
EAN of inner box	4038917905954
EAN case	4038917905961

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

