



### **Product description**

Order number	72.709
Product description	Reaction tube, working volume: 2 ml, material: PP, transparent, without cap, with skirted base, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

### **Size**

Diameter	10.8 mm
Length of product	39 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	210000
Depth of case	388 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0367 cbm
Weight of product	0.0012 kg
Weight of case	6.44 kg
EAN of inner box	4038917083591
EAN case	4038917002493

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

