



Product description

| | |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number | 72.692.405 |
| Product description | Screw cap micro tube, working volume: 1.5 ml, conical base, with knurling, transparent, cap: natural, cap assembled, PCR Performance Tested, 100 piece(s)/bag |

Product characteristics

| | |
|--------------------------|--------------------|
| Graduation | no |
| Cap | standard screw cap |
| Knurls | yes |
| Base shape | conical base |
| Volume of work | 1,5 ml |
| Centrifugation max (RCF) | 20000 x g |

Size

| | |
|-------------------|---------|
| Diameter | 10.8 mm |
| Length of product | 44 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|----------------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | Polypropylene (PP) |
| Colour of cap | natural |
| O-ring material | Ethylene propylene diene rubber (EPDM) |

Purity & certification

| | |
|---------------------------|-------------------------------|
| Satisfies the requirement | IATA |
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Purity standard | PCR Performance Tested |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 110000 |
| Depth of box | 211 mm |
| Width of box | 191 mm |
| Height of box | 141 mm |
| Depth of case | 398 mm |
| Width of case | 223 mm |
| Height of case | 160 mm |
| Case volume | 0.0142 cbm |
| Weight of product | 0.0014 kg |
| Weight of case | 1.94 kg |
| EAN of inner box | 4038917223898 |
| EAN case | 4038917223881 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

