



Product description

| | |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number | 72.730.306 |
| Product description | Screw cap micro tube, working volume: 0.5 ml, skirted conical base, with knurling, transparent, cap: natural, cap enclosed, with printed writing space, sterile, 100 piece(s)/bag |

Product characteristics

| | |
|--------------------------|----------------------|
| Label/ Print | with print |
| Colour of print/label | white |
| Graduation | no |
| Cap | standard screw cap |
| Knurls | yes |
| Base shape | skirted conical base |
| Volume of work | 0,5 ml |
| Centrifugation max (RCF) | 20000 x g |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|---------|
| Diameter | 10.8 mm |
| Length of product | 44 mm |

Material & colours

| | |
|-------------------|----------------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | Polypropylene (PP) |
| Colour of cap | natural |
| O-ring material | Ethylene propylene diene rubber (EPDM) |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | IATA |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 96000 |
| Depth of case | 278 mm |
| Width of case | 348 mm |
| Height of case | 140 mm |
| Case volume | 0.0135 cbm |
| Weight of case | 1.79 kg |
| EAN of inner box | 4038917996556 |
| EAN case | 4038917996563 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

