



Product description

Order number	72.730.711
Product description	Screw cap micro tube, working volume: 0.5 ml, skirted conical base, with knurling, transparent, without cap, sterile, 500 piece(s)/double bag

Product characteristics

Graduation	no
Knurls	yes
Base shape	skirted conical base
Volume of work	0,5 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	10.8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Satisfies the requirement IATA
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 2500
Type of smallest subpackaging double bag
Piece(s) / inner box 500
Piece (s) / outer case 2500
Piece(s) / pallet 130000
Depth of case 598 mm
Width of case 378 mm
Height of case 135 mm
Case volume 0.0305 cbm
Weight of product 0.0012 kg
Weight of case 3.58 kg
EAN of inner box 4038917119566
EAN case 4038917043922

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

