



Product description

Order number	72.733.201
Product description	Screw cap micro tube, for volume limitation to 100 μ l, working volume: 0.5 ml, conical base, with knurling, transparent, without cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Knurls	yes
Example of use	for volume limitation to 100 μ l
Base shape	conical base
Volume of work	0,5 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Satisfies the requirement IATA
Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 500
Type of smallest subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 5000
Piece(s) / pallet 210000
Depth of case 388 mm
Width of case 378 mm
Height of case 250 mm
Case volume 0.0367 cbm
Weight of product 0.0008 kg
Weight of case 4.48 kg
EAN of inner box 4038917366625
EAN case 4038917366618

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

