



Product Specification
**MEDISPONGE, impregnated with 4%
chlorhexidine gluconate solution**

Page 1



Product description

Order number	54.9912.100
Product description	MEDISPONGE, surgical scrub sponge for pre-operative cleaning of hands and fingernails, impregnated with 4% chlorhexidine gluconate solution, material: PUR, (LxWxH): 130 x 63 x 30 mm

Product characteristics

Preparation concentration	4% chlorhexidine gluconate solution
Application	surgical scrub sponge for pre-operative cleaning of hands and fingernails

Size

Width of product	63 mm
Height of product	30 mm
Length of product	130 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-11-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**MEDISPONGE, impregnated with 4%
chlorhexidine gluconate solution**

Page 2

Material & colours

Product material Polyurethane (PUR)
Colour of product white

Purity & certification

Latex-free yes
Pyrogen-free no
Safety data sheet yes

Packaging

Minimum order qty. 450
Type of smallest case
subpackaging
Piece(s) / inner box 75
Piece (s) / outer case 450
Piece(s) / pallet 7200
Depth of case 528 mm
Width of case 465 mm
Height of case 526 mm
Case volume 0.1283 cbm
Weight of case 7.9 kg
EAN of inner box 4038917150491
EAN case 4038917036856

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-11-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT