



### Product description

Order number 87.7709.033  
Product description Mains adapter Mains adapter for AST-2 tube stripper

### Product characteristics

Application Obtaining blood and blood products  
Supply voltage 240 V

### Purity & certification

Product category no medical device | IVD  
CE certified CE supplier  
Batched no

### Packaging

Minimum order qty. 1  
Type of smallest case  
subpackaging  
Piece(s) / inner box 1  
Piece (s) / outer case 1  
Piece(s) / pallet 432  
Depth of case 188 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Mains adapter**

**Page 2**

Width of case	188 mm
Height of case	100 mm
Case volume	0.0035 cbm
Weight of product	0.365 kg
Weight of case	0.54 kg
EAN of inner box	4038917141192
EAN case	4038917067164

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**