



### **Product description**

Order number	97.8800.200
Product description	Stripping / Sealing Tongs, for stripping, segmenting and sealing blood collection tubes

### **Product characteristics**

Application	Obtaining blood and blood products
-------------	------------------------------------

#### **Size**

Width of product	5 cm
Length of product	16 cm

### **Material & colours**

Colour of product	red
-------------------	-----

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	432
Depth of case	188 mm
Width of case	188 mm
Height of case	100 mm
Case volume	0.0035 cbm
Weight of product	0.16 kg
Weight of case	0.26 kg
EAN of inner box	4038917975988
EAN case	4038917975995

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

