



### **Product description**

Order number	74.5340.500
Product description	TUR-BAG, Urine drainage system, 4 l, with pump ball, sterile, tube length: 100 cm, 1 piece(s)/blister

### **Product characteristics**

Colour of print/label	blue
Design	with pump ball
Application	Urine drainage
Volume of work	4 l

### **Size**

Sample volume	4 l
Max. volume	4 l
Diameter of tube	10.1 mm
Length of tube	100 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



### **Material & colours**

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

### **Purity & certification**

Product category	Medical device
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

### **Packaging**

Minimum order qty.	10
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	10
Piece(s) / pallet	560
Depth of case	395 mm
Width of case	295 mm
Height of case	230 mm
Case volume	0.0268 cbm
Weight of product	0.2455 kg
Weight of case	2.46 kg
EAN of inner box	6955824024395
EAN case	6955824024012

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

