



Product description

Order number 74.5685.001
Product description Urine measuring system exchange bag, 3 l, sterile, 1 piece(s)/blister

Product characteristics

Colour of print/label blue
Application Urine drainage
Volume of work 3 l

Size

Sample volume 3 l
Max. volume 3 l

Material & colours

Product material Polyvinylchloride (PVC)
Colour of product transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	Medical device
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	50
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	50
Piece(s) / pallet	3200
Depth of case	395 mm
Width of case	195 mm
Height of case	200 mm
Case volume	0.0154 cbm
Weight of product	0.0523 kg
Weight of case	3.47 kg
EAN of inner box	6955824024456
EAN case	6955824024074

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

