



### Product description

|                     |  |
|---------------------|--|
| Order number        | 74.5228.005  |
| Product description | Urine bag, 2 l, with anti-reflux valve, with drainage valve, sterile, tube length: 90 cm, 1 piece(s)/blister |
| PZN number          | 16333726   |

### Product characteristics

|                       |                     |
|-----------------------|---------------------|
| Colour of print/label | blue                |
| Design                | with drainage valve |
| Application           | Urine drainage      |
| Volume of work        | 2 l                 |
| Anti-reflux valve     | yes                 |

### Size

|                  |        |
|------------------|--------|
| Sample volume    | 2 l    |
| Max. volume      | 2 l    |
| Diameter of tube | 6.6 mm |
| Length of tube   | 90 cm  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine bag, 2 l, with anti-reflux valve, with drainage valve, sterile**

**Page 2**

**Material & colours**

Product material Polyvinylchloride (PVC)  
Colour of product transparent

**Purity & certification**

Product category Medical device  
CE certified CE supplier  
Purity standard sterile  
Sterilisation Gas sterilisation  
Batched yes

**Packaging**

Minimum order qty. 100  
Type of smallest blister  
subpackaging  
Piece(s) / inner box 1  
Piece (s) / outer case 100  
Piece(s) / pallet 4800  
Depth of case 305 mm  
Width of case 390 mm  
Height of case 255 mm  
Case volume 0.0303 cbm  
Weight of product 0.042 kg  
Weight of case 4.5 kg  
EAN of inner box 6971243661041  
EAN case 6971243660808

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-09-01

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**