



### Product description

Order number	74.5282.101
Product description	Secretion bag, 2 l, universal step connector, tube length: 90 cm, 200 piece(s)/case
PZN number	16333583

### Product characteristics

Colour of print/label	blue
Design	universal step connector
Application	Drainage of exudate
Volume of work	2 l
Anti-reflux valve	no

### Size

Sample volume	2 l
Max. volume	2 l
Diameter of tube	6.6 mm
Length of tube	90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



### **Material & colours**

Product material Polyvinylchloride (PVC)  
Colour of product transparent

### **Purity & certification**

Product category Medical device  
CE certified CE supplier  
Batched yes

### **Packaging**

Minimum order qty. 200  
Type of smallest case  
subpackaging  
Piece(s) / inner box 200  
Piece (s) / outer case 200  
Piece(s) / pallet 9600  
Depth of case 580 mm  
Width of case 390 mm  
Height of case 150 mm  
Case volume 0.0339 cbm  
Weight of product 0.035 kg  
Weight of case 7.4 kg  
EAN of inner box 6971243660976  
EAN case 6971243660662

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

