



### **Product description**

Order number	74.5284.101
Product description	Urine bag, 2 l, with anti-reflux valve, tube length: 90 cm, 200 piece(s)/case
PZN number	16333749

### **Product characteristics**

Colour of print/label	blue
Application	Urine drainage
Volume of work	2 l
Anti-reflux valve	yes

### **Size**

Sample volume	2 l
Max. volume	2 l
Diameter of tube	6.6 mm
Length of tube	90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polyvinylchloride (PVC)  
Colour of product transparent

### **Purity & certification**

Product category Medical device  
CE certified CE supplier  
Batched yes

### **Packaging**

Minimum order qty. 200  
Type of smallest case  
subpackaging  
Piece(s) / inner box 200  
Piece (s) / outer case 200  
Piece(s) / pallet 9600  
Depth of case 580 mm  
Width of case 390 mm  
Height of case 150 mm  
Case volume 0.0339 cbm  
Weight of product 0.035 kg  
Weight of case 7.4 kg  
EAN of inner box 6971243660983  
EAN case 6971243660679

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

