



### Product description

Order number 79.8710.575  
Product description Paper roll, Consumables for protocol printer, white

### Product characteristics

Application Warming

### Size

Width of product 7.5 cm  
Length of product 8.3 cm

### Material & colours

Product material paper  
Colour of product white

### Purity & certification

Product category no medical device | IVD  
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Packaging**

Minimum order qty.	10
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	50
Piece(s) / pallet	3000
Depth of case	400 mm
Width of case	400 mm
Height of case	170 mm
Case volume	0.0272 cbm
Weight of product	0.269 kg
Weight of case	14.15 kg
EAN of inner box	4038917975865
EAN case	4038917975872

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

