



### **Product description**

Order number	79.8710.610
Product description	Adaptation compress TSC, Replacement gel compress for thawing out cryopreserved blood stem cell preparations in the SAHARA TSC

### **Product characteristics**

Application	Warming
-------------	---------

### **Size**

Width of product	16 cm
Length of product	23.5 cm

### **Purity & certification**

Product category	no medical device   IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Packaging**

Minimum order qty.	5
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	5
Piece(s) / pallet	960
Depth of case	320 mm
Width of case	220 mm
Height of case	140 mm
Case volume	0.0099 cbm
Weight of product	0.488 kg
Weight of case	2.6 kg
EAN of inner box	4038917978545
EAN case	4038917978552

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

