



Product Specification
Ventilator for Blood Gas Monovette® ,
transparent/orange, length: 24.5 mm

Page 1



Product description

Order number 14.1148
Product description Ventilator for Blood Gas Monovette® , transparent/orange,
length: 24.5 mm, accessories for the Blood Gas
Monovette® , 500 piece(s)/case

Product characteristics

Application accessories for the Blood Gas Monovette®
Example of use Ventilation of the Blood Gas -Monovette®
Suitable for Blood Gas Monovette®

Size

Length of product 24.5 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent/orange

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Ventilator for Blood Gas Monovette® ,
transparent/orange, length: 24.5 mm

Page 2

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Latex-free	yes
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	168000
Depth of box	191 mm
Width of box	129 mm
Height of box	129 mm
Depth of case	288 mm
Width of case	218 mm
Height of case	150 mm
Case volume	0.0094 cbm
Weight of product	0.0009 kg
Weight of case	1.34 kg
EAN of inner box	4038917085465
EAN case	4038917075312

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT