



### Product description

Order number	19.930.125
Product description	Blood gas capillaries, preparation: calcium-balanced lithium heparin, nominal volume: 125 µl, (LxØ): 110 x 2.05 mm, 250 piece(s)/container

### Product characteristics

Type of collection	capillary
Type of preparation	calcium-balanced lithium heparin
Preparation concentration	70 IU/ml blood
Application	blood gas analysis

### Size

Sample volume	125 µl
Diameter	2.05 mm
Length of product	110 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polyethylene terephthalate (PET)  
Colour of product transparent

### **Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Latex-free yes  
Batched yes

### **Packaging**

Minimum order qty. 250  
Type of smallest container  
subpackaging  
Piece(s) / inner box 250  
Piece (s) / outer case 16000  
Piece(s) / pallet 1280000  
Depth of case 378 mm  
Width of case 288 mm  
Height of case 180 mm  
Case volume 0.0196 cbm  
Weight of product 0.0003 kg  
Weight of case 6.64 kg  
EAN of inner box 4038917234283  
EAN case 4038917234276

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

