



### **Product description**

Order number	19.930.175
Product description	Blood gas capillaries, preparation: calcium-balanced lithium heparin, nominal volume: 175 µl, (LxØ): 100 x 2.3 mm, 200 piece(s)/container

### **Product characteristics**

Type of collection	capillary
Type of preparation	calcium-balanced lithium heparin
Preparation concentration	70 IU/ml blood
Application	blood gas analysis

### **Size**

Sample volume	175 µl
Diameter	2.3 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polyethylene terephthalate (PET)  
Colour of product transparent

### **Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Latex-free yes  
Batched yes

### **Packaging**

Minimum order qty. 200  
Type of smallest container  
subpackaging  
Piece(s) / inner box 200  
Piece (s) / outer case 12800  
Piece(s) / pallet 921600  
Depth of case 333 mm  
Width of case 303 mm  
Height of case 190 mm  
Case volume 0.0192 cbm  
Weight of product 0.0003 kg  
Weight of case 5.5 kg  
EAN of inner box 4038917234344  
EAN case 4038917234337

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

