



Product description

Order number	51.931.175
Product description	Blood gas capillary set, preparation: calcium-balanced lithium heparin, nominal volume: 175 µl, (LxØ): 100 x 2.3 mm, set components: 800 PET plastic capillaries, 175 µl (4 x 200 pieces REF 19.930.175), 1,600 blue quick-release caps (4 x 400 pieces REF 65.935.230) and 1,000 mixing wires (4 x 250 pieces REF 95.936), 1 piece(s)/case

Product characteristics

Type of collection	capillary
Type of preparation	calcium-balanced lithium heparin
Preparation concentration	70 IU/ml blood
Application	blood gas analysis

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Blood gas capillary set, calcium-balanced lithium
heparin, 175 µl**

Page 2

Size

Sample volume	175 µl
Diameter	2.3 mm
Length of product	100 mm

Material & colours

Product material	Polyethylene therephthalate (PET)
Closure material	Thermoplastic elastomers (TPE)

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	16
Piece(s) / pallet	432
Depth of case	778 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0558 cbm
Weight of product	0.6551 kg
Weight of case	14.55 kg
EAN of inner box	4038917234467
EAN case	4038917234450

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT