



Product description

| | |
|---------------------|---|
| Order number | 95.936 |
| Product description | Mixing wire for blood gas capillaries, (LxØ): 9 x 1 mm, accessories for the Blood Gas Monovette® , 250 piece(s)/container |

Product characteristics

| | |
|--------------------|--|
| Type of collection | capillary |
| Application | accessories for the Blood Gas Monovette® |
| Suitable for | blood gas capillaries |

Size

| | |
|-------------------|------|
| Diameter | 1 mm |
| Length of product | 9 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Mixing wire for blood gas capillaries, (LxØ): 9 x 1 mm

Page 2

Material & colours

Product material metal

Purity & certification

Product category In vitrodiagnostic device, CE

CE certified CE

Batched yes

Packaging

Minimum order qty. 250

Type of smallest container

subpackaging

Piece(s) / inner box 250

Piece (s) / outer case 41250

Piece(s) / pallet 412500

Depth of case 268 mm

Width of case 191 mm

Height of case 100 mm

Case volume 0.0051 cbm

Weight of product 0.0001 kg

Weight of case 2.43 kg

EAN of inner box 4038917234610

EAN case 4038917234603

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT