



Product Specification
S-Sedivette® , 3.5 ml, cap black, (LxØ): 130 x 8 mm, with paper label

Page 1



Product description

Order number	06.1690.100
Product description	S-Sedivette® , blood sedimentation, 4NC, preparation: Trisodium citrate/citric acid buffer solution, nominal volume: 3.5 ml, (LxØ) with cap: 146 x 8 mm, with paper label, label/print: black, cap black, colour code ISO, sterile, 50 piece(s)/case

Product characteristics

Type of collection	venous
Type of preparation	Trisodium citrate/citric acid buffer solution
Preparation concentration	0.105 mol/l pH 5.5
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	black
Graduation	no
Cap	membrane screw cap
Closure type	screw cap
Application	blood sedimentation

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-11-15

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Size

Sample volume	3.5 ml
Diameter	8 mm
Length including cap	146 mm
Length excluding cap	130 mm
Length including cap and plunger	167 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	black
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	36000
Depth of box	143 mm
Width of box	171 mm
Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Sedivette® , 3.5 ml, cap black, (LxØ): 130 x 8 mm, with paper label

Page 3

Height of case	190 mm
Case volume	0.0192 cbm
Weight of product	0.0052 kg
Weight of case	3.29 kg
EAN of inner box	4038917133128
EAN case	4038917024716

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-11-15

This document was prepared by EDP support and is valid without signature.



SARSTEDT