



### **Product description**

Order number	85.1026
Product description	Incision safety lancet, Safety-Heel® Newborn, light blue/white, penetration depth: 1 mm, blade width: 2.5 mm, for heel puncture, 50 piece(s)/case

### **Product characteristics**

Example of use	Blood collection from newborn babies
Puncture site	Heel
Penetration depth	1 mm
Safety product	yes

### **Size**

Width of product	35 mm
Length of product	40 mm
Width of blade	2.5 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Colour of product light blue/white

### **Purity & certification**

Product category Medical device  
CE certified CE supplier  
Purity standard sterile  
Sterilisation Gamma irradiation  
Batched yes

### **Packaging**

Minimum order qty. 200  
Type of smallest case  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 200  
Piece(s) / pallet 96000  
Depth of case 273 mm  
Width of case 262 mm  
Height of case 160 mm  
Case volume 0.0114 cbm  
Weight of product 0.0094 kg  
Weight of case 2.08 kg  
EAN of inner box 6945630110474  
EAN case 6945630110535

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

