



### Product description

|                     |   |
|---------------------|---|
| Order number        | 15.1675.100   |
| Product description | Multivette® 600 Lithium heparin gel LH, venous and capillary blood samples, preparation: Lithium heparin gel, nominal volume: 600 µl, (LxØ) with cap: 53.3 x 8 mm, with print, cap green, screw cap, flat base, 100 piece(s)/case |

### Product characteristics

|                            |                                    |
|----------------------------|------------------------------------|
| Type of collection         | Venous and capillary               |
| Type of preparation        | Lithium heparin gel                |
| Colour code                | ISO                                |
| Label/ Print               | with print                         |
| Colour of print/label      | green                              |
| Closure type               | screw cap                          |
| Application                | venous and capillary blood samples |
| Base shape                 | flat base                          |
| Recommended centrifugation | 10000 x g                          |
| Centrifugation time        | 5 Min                              |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### Size

|                      |         |
|----------------------|---------|
| Sample volume        | 600 µl  |
| Diameter             | 8 mm    |
| Length including cap | 53.3 mm |
| Length excluding cap | 49.8 mm |

### Material & colours

|                   |                                   |
|-------------------|-----------------------------------|
| Product material  | Polypropylene (PP)                |
| Colour of product | transparent                       |
| Closure material  | High Density Polyethylene (HD-PE) |
| Colour of cap     | green                             |

### Purity & certification

|                  |                                |
|------------------|--------------------------------|
| Product category | In vitro diagnostic device, CE |
| CE certified     | CE                             |
| Batched          | yes                            |

### Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1000          |
| Type of smallest subpackaging | case          |
| Piece(s) / inner box          | 100           |
| Piece (s) / outer case        | 1000          |
| Piece(s) / pallet             | 63000         |
| Depth of box                  | 201 mm        |
| Width of box                  | 81 mm         |
| Height of box                 | 81 mm         |
| Depth of case                 | 423 mm        |
| Width of case                 | 213 mm        |
| Height of case                | 185 mm        |
| Case volume                   | 0.0167 cbm    |
| Weight of product             | 0.0024 kg     |
| Weight of case                | 3.01 kg       |
| EAN of inner box              | 4038917133005 |
| EAN case                      | 4038917024648 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

