



**Product Specification**  
**Microvette® CB 300 Serum CAT, 300 µl, cap white, push cap, flat base**

**Page 1**



### Product description

Order number	16.440
Product description	Microvette® CB 300 Serum CAT, Veterinary medicine, preparation: Clotting activator, nominal volume: 300 µl, (LxØ) with cap: 46.5 x 10.8 mm, cap white, push cap, flat base, inner receptacle: trumpet-shaped, 100 piece(s)/bag

### Product characteristics

Type of collection	capillary
Type of preparation	Clotting activator
Colour code	EU
Cap	Soft stopper with retaining loop
Closure type	push cap
Application	Veterinary medicine
Base shape	flat base
Recommended centrifugation	10000 x g
Centrifugation time	5 Min
Shape of inner tube	trumpet-shaped

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Microvette® CB 300 Serum CAT, 300 µl, cap**  
**white, push cap, flat base**

**Page 2**

**Size**

Sample volume	300 µl
Diameter	10.8 mm
Length including cap	46.5 mm
Length excluding cap	43.7 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	white

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	80000
Depth of case	378 mm
Width of case	288 mm
Height of case	180 mm
Case volume	0.0196 cbm
Weight of product	0.002 kg
Weight of case	2.4 kg
EAN of inner box	4038917268271
EAN case	4038917282703

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**