



## Product Specification

**Minivette® POCT Lithium heparin LH, 50 µl,  
plunger orange, colour code EU, 200 piece(s)/bag**

**Page 1**



### Product description

Order number	17.2112.050
Product description	Minivette® POCT Lithium heparin LH, nominal volume: 50 µl, preparation: Heparin, orange, colour code EU, 200 piece(s)/bag, 200 piece(s)/case

### Product characteristics

Type of collection	capillary
Type of preparation	Heparin
Colour code	EU
Application	capillary blood collection
Coding colour	orange

### Size

Sample volume	50 µl
---------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Minivette® POCT Lithium heparin LH, 50 µl,  
plunger orange, colour code EU, 200 piece(s)/bag**

**Page 2**

### Material & colours

Product material	Cellulose acetate (CA)
Colour of product	transparent

### Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Latex-free	yes
Batched	yes

### Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	200
Piece (s) / outer case	2000
Piece(s) / pallet	144000
Depth of box	140 mm
Width of box	99 mm
Height of box	117 mm
Depth of case	513 mm
Width of case	293 mm
Height of case	140 mm
Case volume	0.021 cbm
Weight of product	0.001 kg
Weight of case	2.68 kg
EAN of inner box	4038917310345
EAN case	4038917310338

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**