



Product description

| | |
|---------------------|---|
| Order number | 20.2331.100 |
| Product description | Microvette® APT 250 EDTA K2E, capillary blood collection, preparation: K2 EDTA, K2E, nominal volume: 250 µl, Ø: 13 mm, cap violet, cap, with pierceable membrane, round base, inner receptacle: cylindrical, incl. quick-release cap, non-sterile, 50 piece(s)/case |

Product characteristics

| | |
|-----------------------|----------------------------|
| Type of collection | capillary |
| Type of preparation | K2 EDTA |
| Colour code | ISO |
| Colour of print/label | transparent/violet |
| Cap | membrane cap |
| Closure type | cap |
| Application | capillary blood collection |
| Base shape | round base |
| Volume of work | 250 µl - 250 µl |
| Shape of inner tube | cylindrical |

Size

| | |
|----------------------|--------|
| Sample volume | 250 µl |
| Diameter | 13 mm |
| Length excluding cap | 75 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | Polyethylene (PE) |
| Colour of cap | violet |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | P 650, ADR |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Purity standard | non-sterile |
| Batched | yes |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Microvette® APT 250 EDTA K2E, 250 µl, cap
violet, cap, round base**

Page 2

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 32000 |
| Depth of box | 201 mm |
| Width of box | 128 mm |
| Height of box | 72 mm |
| Depth of case | 378 mm |
| Width of case | 268 mm |
| Height of case | 221 mm |
| Case volume | 0.0224 cbm |
| Weight of product | 0.0075 kg |
| Weight of case | 4.01 kg |
| EAN of inner box | 4038917946070 |
| EAN case | 4038917946087 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT