



Product description

Order number 95.1048
Product description Capillary holder, capillary blood collection, 5 piece(s)/bag

Product characteristics

Type of collection capillary
Application capillary blood collection

Material & colours

Colour of product red

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 5
Type of smallest subpackaging bag
Piece(s) / inner box 5
Piece (s) / outer case 1000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Capillary holder

Piece(s) / pallet	108000
Depth of case	308 mm
Width of case	223 mm
Height of case	195 mm
Case volume	0.0134 cbm
Weight of product	0.0019 kg
Weight of case	2.31 kg
EAN of inner box	4038917085120
EAN case	4038917003810

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

